


Notes:  
 NC channels labeled as 'X' are reserved for future applications including extra channels in 192-channel probes.  
 Ain1 to Ain128 are transducer channels 1 to 128.  
 Use Ain1 to Ain64 for 64-channel probes. Use Ain1 to Ain96 for 96-channel probes.



J1 Cannon QLC, (part number Cannon 127050-0530, the same as TE 6565204) is the probe input connector.

Title		
TE 6565204 / 6565199 / Cannon QLC260 PhotoSound MoleculUS128		
Size	Document Number	Rev
B	This is a standard PhotoSound probe pinout starting 2021	1.0
Date:	Tuesday, January 12, 2021	Sheet 1 of 3